

RAPID DETERMINATION OF ACTINIDES IN WATER SAMPLES USING TEVA AND TK221 RESIN CARTRIDGES

*Nóra Vajda – István Papp**
RADANAL Ltd., Budapest, Hungary
*ISOTOPTECH Co., Debrecen, Hungary

in co-operation with

Steffen Happel
TRISKEM Int. SAS, Bruz, France

eUGM, Nov 24, 2021

Contents

- **Properties of TEVA and TK221 resins**
- **Models:**
 - Separation of actinides on TEVA**
 - Separation of actinides on TK221**
- **Determination of actinides in tap water**
 - Determination of actinides in sea water**

Goal

to develop combined procedure for the selective separation of Pu, Am-Cm, U, Np, Th using **chromatographic resins**.

The procedure should be

- adequate for water sample
- adequate for measurement by α spectrometry
- rapid
- robust

to analyse actinides of elevated activities

Analytical applications:

process water

contaminated water (drinking, surface, ground, sea) in emergency situation.

Goal

to develop combined procedure for the selective separation of Pu, Am-Cm, U, Np, Th using **chromatographic resins**.

The procedure should be

- adequate for water sample
- adequate for measurement by α spectrometry
- rapid
- robust

to analyse actinides of elevated activities

Analytical applications:

process water

contaminated water (drinking, surface, ground, sea) in emergency situation.

Resins of high distribution ratios and selectivities for actinides are needed.

Properties of actinides on specific resins

Single column option?

one resin retains all actinides

separation is

sensitive to the conditions

Multi-column option?




resins retain certain actinides

more robust

k' capacity factors of An's on EC resins

in HNO_3


Pu
Np
Th
Am
U




k' capacity factors of An's on EC resins

in 3M HNO₃


TEVA



TRU





n-DGA



k' capacity factors of An's on TEVA resin

in HNO_3 and HCl



Load

from 3M HNO_3 :

Pu(IV), Np(IV), Th

→ effluent: U(VI), Am(III)

Elution possibilities:

from 9M HCl :

Th

from 0.1M HCl :

Pu(IV), Np(IV)

from 9M $\text{HCl}/\text{Ti}^{3+}$:


Pu(III)

k' capacity factors of An's on TRU resin


in HNO_3 and HCl

~~Pu~~
~~Np~~
~~Th~~
 (Am)
 U

TRU from 3M HNO_3



TRU from HCl



After removal of Pu, Np, Th

Load

from 3M HNO_3 :

U(VI), (Am)

→ effluent: (Am)

from 0.25M HCl :

Am

Elution possibilities:

from 0.1M $\text{NH}_4\text{HC}_2\text{O}_4$:


U(VI)

k' capacity factors of An's on DGA resin


from HNO_3 and HCl

~~Pu~~
~~Np~~
~~Th~~
 Am
~~U~~

DGA from HNO_3



DGA from HCl



After removal of Pu, Np, Th, U

Load

from 3M HNO_3 :
 Am(III), (U)
 → effluent: (U)

Elution possibilities:

from 0.25M HCl :
 Am

TK221 resin

Combination of TRU resin and DGA resin

Ln's

Designed for

- separation of **lanthanides**, e.g. ^{177}Lu
- separation of **actinides**

Pu

Np

Th

Am

U

Out of natural matrix elements (alkali, alkaline earth, Al) only Ca is slightly retained from HNO_3 .

Many III and IV valent elements (e.g. Fe, Zr) are well retained from HNO_3 .

Product sheet: TK221 Resin, https://www.triskem-international.com/scripts/files/60cb473fbbe7e8.32187398/PS_TK221-Resin_EN_210607.pdf

S. Happel: "An overview over some new extraction chromatographic resins and their application in radiopharmacy" presented on the 4th of June 2019 at the 102nd Canadian Chemistry Conference and Exhibition (CCCE 2019) in Quebec City, QC


N. Vajda et al. (RadAnal): "Report: Investigation on TEVA/TK221 resins for separation of actinides", Budapest, April 2021

Basic concept

Actinides pre-concentration with Ca phosphate co-precipitation

Load:

15 mL 3M HNO₃/
1M Al(NO₃)₃
Oxidation state
adjustment



Rinse:

18 mL 3M HNO₃

Rinse:

10 mL 3M HNO₃

Th strip:

15 mL 9M HCl

Pu-Np strip:

15 mL 0.1M HCl/0.05M
HF/0.03M Ti³⁺



Th/NdF₃

Pu-Np/NdF₃

Rinse:

10 mL 3M HNO₃

Am strip:

15 mL 0.25M HCl

U strip:

15 mL 0.1M NH₄HC₂O₄



Am/NdF₃


U/NdF₃

Basic concept – extended for Pu/Np separation

Actinides pre-concentration with Ca phosphate co-precipitation

Load:

15 mL 3M HNO₃/
1M Al(NO₃)₃
Oxidation state
adjustment



Rinse:

18 mL 3M HNO₃

Rinse:

10 mL 3M HNO₃

Th strip:

15 mL 9M HCl

Pu strip:

15 mL 9M HCl/0.03M Ti³⁺

Np strip:

20 mL 0.1M HCl/
0.1M HF

Rinse:

10 mL 4M HCl/0.1M HF

Am strip:

15 mL 0.25M HCl

U strip:

15 mL 0.1M NH₄HC₂O₄



Th/NdF₃

Pu/NdF₃

Np/NdF₃

Am/NdF₃

U/NdF₃

Optimization of the EC separation

Optimization of Am/U separation on TK221

Model

simulating water sample

Ca phosphate co-precip.
Fe and Ca in the load

Oxidation state adjustment:
sulfamic acid/ascorbic acid/
NaNO₂

Tracers:

individual: ²⁴¹Am, ²³³U
about 50 Bq in each test



Measurement by LSC

Effect of
Ca in load

No Ca



50 mg Ca



Am
is hold back


Optimization of Am/U separation on TK221

Model

simulating water sample

Tracers:

individual: ^{241}Am , ^{233}U




Good Am/U separation by increased eluent volume!

Separation of Am, U, Th, Pu on TEVA/TK221 and TK221

Model
simulating water
sample



Tracers:
individual:
 ^{241}Am , ^{233}U
 ^{230}Th , ^{239}Pu



Separation of Th, Pu(Np), Am, U on TEVA/TK221 and TEVA

Model

simulating water sample, Tracers: individual: ^{241}Am , ^{233}U , ^{230}Th , ^{239}Pu



Results of optimization

Model

simulating water sample

Ca phosphate co-precip.
Fe and Ca in the load

Oxidation state adjustment:
sulfamic acid/ascorbic acid/
NaNO₂

Tracers:

individual: ²³⁰Th, ²³⁹Pu,
²⁴¹Am, ²³³U

about 50 Bq in each test

Measurement by LSC

Recovery

		U %	Am %	Th %	Pu %
TEVA/TK221 cartridges					
load	15 mL load				
wash	18 mL 3M HNO ₃	2-5	0	0	0
TEVA cartridge					
wash	10 mL 3M HNO ₃	0,1-0,2	0	0.1	0
Th strip	15 mL 9M HCl	0.3	0	97	0
Pu strip	15 mL 0.1M HCl/0.05M HF/0.03M TiCl ₃	0.2	0	0	106
TK221 cartridge					
wash	15 mL 4M HCl/0.2M HF	0.5	0.7	0	0
Am strip	30 mL 0.25M HCl	0.8	102	0	0
U strip	15 mL 0.1M NH ₄ HC ₂ O ₄	97	0.7	0	0
	SUM %	101-104	103	97	106

High recoveries, some U contamination (<1%) is possible!

Determination of actinides in water samples

Determination of actinides in water samples

800 mL water:
 TAP or SEA
 spiked with mixed
 tracers:
²⁴¹Am, ²³³U, ²³⁰Th,
²³⁹Pu, ²³⁷Np
 of 0.1-0.5 Bq
 Ca phosphate
 co-precipitation,
 TEVA/TK221
 separation,
 μ-co-precipitation
 α spectrometry

	Actinides determination			
	without Np		with Np separation	
	yield %	yield unc %	yield %	yield unc %
TAP water				
²³⁰ Th	90	8	86	7
²³⁹ Pu	108	7	95	7
²³⁷ Np	-		91	9
²⁴¹ Am	103	7	97	6
²³³ U	103	7	70	7
SEA water				
²³⁰ Th	71	7	61	6
²³⁹ Pu	91	7	87	6
²³⁷ Np	-		93	8
²⁴¹ Am	89	7	92	6
²³³ U	88	7	59	6

Results

Recoveries of all An's are
 acceptable high:
 > 60%

Recoveries in TAP are higher
 than in SEA.

Detectable contamination in
 α sources:
 < 1% ²³⁹Pu in ²³⁷Np source

To get high Pu recovery in
 α source the Pu strip solution
 (9M HCl) has to be evaporated,
 diluted.

Summary

A new procedure for the separation of U, Pu, Am-Cm and Np in water samples has been developed where Actinides are

- *pre-concentrated with Ca phosphate precipitate*
- *chromatographically separated on TEVA/TK221 cartridges,*
- *determined by alpha spectrometry.*

For all actinides

- *chemical recoveries are acceptable high ($\geq 60\%$),*
- *decontamination factors are acceptable high (>100),*
- *resolution of the alpha sources are acceptable good (< 40 keV),*
- *sensitivities are acceptable high (<0.1 Bq/kg)*

*The whole analysis can be done in **1 day**.*