



# TrisKem International

Application of extraction chromatography in the purification and QC of radiometals for use in Nuclear Medicine and Radiopharmacy

S. Happel, M. Bas, I. Dovhyi

07/06/2023

CSC  VANCOUVER  
2023 JUNE 4-8 | #CSC2023



# TrisKem International



- Based in Rennes (France)
- Independent company since 02/07
  - Formerly part of Eichrom Europe
  - ISO 9001 since 2007
- Main product line: extraction chromatographic resins
- Staff : 20
- R&D, QC and TechSupport group:
  - 4 RadChem PhD, 3 Technicians
- R&D: Development of new resins, techniques and applications
- Several domains


Radiopharmacy  
and  
Nuclear Medicine

Environment and  
Bioassay

Geochemistry  
and  
Metals Separation



Decommissioning

- Purification of radiometals for use in diagnostics and therapy
  - Mainly: separation of radionuclides from irradiated targets
    - Diagnostics: Zr-89, Cu-64, Ga-68, Ge-68, Ti-44/5, Tc-99m, Sc-43/4...
    - Therapy: Ac-225, Pb-212, Lu-177, Tb-161, Cu-67, Sn-117m, Sc-47...
  - Also: QC, valorization/recycling and waste/effluent treatment
  - Resin characterization and method development done 'cold'
    - Cooperation with cyclotrons & reactors, universities,...
  - Choice of right resin(s):
    - Rapid separations, high purity of products
    - Radiolysis stability
    - No selectivity for target material, high selectivity for product
    - Elution under 'soft' conditions in small volume => labelling/injection
    - Facile automatization e.g. cassettes / no evaporation steps
  - Combining several resins can facilitate the separation
    - Conversion (high acid to dilute acid)
    - Removal of impurities upfront



# Ga-68 separation from Zn targets

- Irradiation of Zn-68 targets in cyclotron
- Ga-68 separation on ZR Resin
  - No selectivity for Zn (target material)
  - Loading possible from:
    - dilute acid (**liquid targets => typically HNO<sub>3</sub>**)
    - >6M HCl (**solid targets**)
  - Rinse under loading condition
  - Elution with ~1 - 2M HCl
  - Too acidic for injection or labelling



- Ga-68 'conversion' necessary
  - Evaporation & dissolution difficult to automatize
- Easier => use of another resin
- TK200 Resin (TOPO) load from 1 - 2M HCl
- Rinse with e.g. 1 - 2M HCl
- Elution in 2 – 3 BV water, dilute acid,..

⇒ **New IAEA TechDoc:**

<https://www-pub.iaea.org/books/IAEABooks/13484/Gallium-68-Cyclotron-Production>

# Cyclotron production of Ga-68

Rodnick et al. *EJNMMI Radiopharmacy and Chemistry* (2020) 5:25  
<https://doi.org/10.1186/s41181-020-00106-9>


EJNMMI Radiopharmacy  
and Chemistry


RESEARCH ARTICLE

Open Access

Cyclotron-based production of  $^{68}\text{Ga}$ ,  $[^{68}\text{Ga}]\text{GaCl}_3$ , and  $[^{68}\text{Ga}]\text{Ga-PSMA-11}$  from a liquid target



Melissa E. Rodnick<sup>1</sup>, Carina Sollert<sup>2</sup>, Daniela Stark<sup>3</sup>, Mara Clark<sup>1</sup>, Andrew Katsifis<sup>3</sup>, Brian G. Hockley<sup>1</sup>, D. Christian Parr<sup>2</sup>, Jens Frigell<sup>2</sup>, Bradford D. Henderson<sup>1</sup>, Monica Abghari-Gerst<sup>1</sup>, Morand R. Piert<sup>1</sup>, Michael J. Fulham<sup>4</sup>, Stefan Eberl<sup>5</sup>, Katherine Gagnon<sup>2</sup> and Peter J. H. Scott<sup>1</sup> 



**Table 1** High level schemes of  $[^{68}\text{Ga}]\text{GaCl}_3$  purifications

	Scheme A*	Scheme B
① ZR Load	< 0.1 M $\text{HNO}_3$	
② ZR Wash	15 mL 0.1 M $\text{HNO}_3$	
③ ZR Elution / Trapping on TK200	5–6 mL ~ 1.75 M HCl	
④ TK Wash	–	3.5 mL 2.0 M NaCl in 0.13 M HCl
⑤ TK Elution	$\text{H}_2\text{O}$	1–2 mL $\text{H}_2\text{O}$ followed by dilute HCl to formulate

\*Process as reported previously (Nair et al. 2017)

## • Solid targets:

- J. Kumlin et al.
- ZR, LN & TK200 for solid targets
  - High Ga-68 activities
  - ARTMS/Odense: 10 Ci production
- One column separation possible using TK400 Resin => solid targets
  - Ga retention on TK400 from high HCl
  - No Zn retention
  - Faster flow than ZR Resin
- W. Tieu et al. use of single TK400 cartridge for solid Zn targets
- Svedjehed et al. use of TK400/A8/TK200 for solid Zn targets



Demystifying solid targets: Simple and rapid distribution-scale production of  $[^{68}\text{Ga}]\text{GaCl}_3$  and  $[^{68}\text{Ga}]\text{Ga-PSMA-11}$

Johan Svedjehed, Martin Pärnaste, Katherine Gagnon\*

Cyclotrons and TRACERcenter, GEMS PET Systems AB, GE Healthcare, Uppsala, Sweden

# TK400 Resin

- Long chained alcohol – initial work by A. Knight et al.
- Retention only at high HCl concentration, elution in low HCl, water,...
- Main application: **Pa separation**
- Other applications:
  - Also retains Mo, Nb, **Fe, Ga, Po**
  - **Fe separation** (higher cap than TRU)
  - **Ga removal from Cu-67 & Fe, Nb removal from Zr**
  - **Ga separation from Zn at high HCl**
  - Nb separation from Zr possible (Nb-90)
  - Under further testing for At separation (elution...)
  - Bleeding => upcoming TK401 Resin





- Preferably Ga elution in 1.5 – 2M HCl => Fe remains retained on TK400
- Ga elution onto TK200

# Cu-61/4 separation on TK201

## ➤ Cu-61/4 separation from solid Ni targets

- Ni targets dissolved in high HCl
    - CU Resin generally not used
  - Typically separation via anion exchange
    - Elevated elution volumes
  - Use of TK201 instead (sharper elution)
    - Originally for Tc separation
    - No selectivity for Ni, good Cu retention, Zn very well retained
  - Load and rinse at 6M HCl
    - Ni removal and recovery/recycling
  - Co elution with 4 – 5M HCl
  - First tests: Cu elution with 0.5M HCl
    - Zn remains retained
    - Ga and Fe separation
- ⇒ Eluate too acidic - requires further treatment => requires improvement



# Cu-61/4 separation on TK201

## ➤ Improvement:

- Use of TBP (or TK400) upfront for Fe/Ga removal
  - allows for Cu elution in 0.05M HCl => suitable for labeling
  - Zn remains retained
- Problem: TK201 'acid sponge' => eluate higher than 0.05M HCl
- Gagnon et al. use of NaCl/HCl for better pH control of eluate

- TK201 also used with CU Resin (and TK400) for Cu-67 separation from Zn targets
  - Conversion from 6M HCl to 0.05M HCl and Zn removal
- Currently being tested for Co separation (2x TK201) from Ni targets and for Zn separation

Svedjehed et al. *EJNMMI Radiopharmacy and Chemistry* (2020) 5:21  
<https://doi.org/10.1186/s41181-020-00108-7>


EJNMMI Radiopharmacy  
and Chemistry

RESEARCH ARTICLE

Open Access

## Automated, cassette-based isolation and formulation of high-purity [<sup>61</sup>Cu]CuCl<sub>2</sub> from solid Ni targets


Johan Svedjehed<sup>1</sup>, Christopher J. Kutryreff<sup>2</sup>, Jonathan W. Engle<sup>2,3</sup> and Katherine Gagnon<sup>1\*</sup>





# Tb-161 separation

- nca Lu-177 more frequently used but strong interest in nca Tb-161
- Part of the ‘Swiss knife of nuclear medicine’ => Tb isotopes



Tb 149		Tb 152		Tb 155	Tb 161
4.2m	4.1h	4.2m	17.5h	5.32d	6.90d
ε	ε	γ283;	ε	ε	β <sup>-</sup> 0.5; 0.6...
β <sup>+</sup>	α3.97	160...	β <sup>+</sup> 2.8...	γ87;	γ 26; 49; 75...
α3.99	β <sup>+</sup> 1.8	ε; β <sup>+</sup> ...	γ 344;	105...	e <sup>-</sup>
γ796;	γ352;	γ344;	586;	180, 262	
165...	165...	411...	271...		

## Terbium: a new ‘Swiss army knife’ for nuclear medicine

Source: <https://cerncourier.com/a/terbium-a-new-swiss-army-knife-for-nuclear-medicine/>



- Typically irradiation of several hundreds of mg (or more)
- Work on upscale on-going (incl. recycling and decontamination of effluents)  
=> originally 0.5 - 1g targets / now: >1g for Tb-161 and >5 – 10g for Lu-177
- Separation via TK211/2/3 Resins
  - Mixture of extractants (HDEHP, HEH[EHP], Cyanex 272, Cyanex 572)
  - Inert support containing aromatic groups + higher capacity for extractant
  - Long-chained alcohol as radical scavenger
  - 30μm particles => smaller (10 – 15μm) possible but high pressure drop...<sub>9</sub>

# Tb separation from 1000 mg Gd targets

- 500 – 1000 mg Gd in cold work, hot 200 – 300 mg Gd targets
- Irradiated target typically dissolved in >3M HNO<sub>3</sub> => Conversion via TK221 Resin
- **Sequential separation** on TK212/TK211
- Final conversion to dilute HCl on TK221 + trace nitrate removal on AIX
- Mainly Tb-161, also Tb-155


A simple and automated method for <sup>161</sup>Tb purification and ICP-MS analysis of <sup>161</sup>Tb

Scott W. McNeil<sup>1</sup>, Michiel Van de Voorde<sup>2</sup>, Chengcheng Zhang<sup>3</sup>, Maarten Ooms<sup>2</sup>, François Bénard<sup>3</sup>, Valery Radchenko<sup>1,4</sup> and Hua Yang<sup>1,5\*</sup>




# Tb separation from 1000 mg Gd targets

- Initial separation on TK212 – 147 mL column (30cm x 2.5cm)
- Flow rate ~15 mL/min
- Gd recovery => very expensive & difficult to find
- Tb separation from Gd and Dy – ideally using online detection
- Fine purification on TK211 (29 mL)



Tb separation from 1000 mg Gd on TK212 (147 mL column)

# Tb purification on TK211




- Direct load of Tb fraction from TK212 onto TK211 (29 mL – 30cm x 1.1cm)
- Flow rate ~15 mL/min
- Gd breakthrough during load & rinse with 0.5M HNO<sub>3</sub> (alternatively HCl)
- Tb elution (Dy sufficiently well removed before) preferably in **>3M HNO<sub>3</sub>**
- Conversion to dilute HCl via TK221, A8 for nitrate removal
- **Option: LN3 cartridge for Dy removal from Tb before use (e.g. after long shipment)**

- Ac-225 separation chemistry well established
- Currently typically DGA (mainly B) is used fore Ac/Ra separation
  - Open questions:
    - Imperfect La/Ac separation (suggested: additional separation on LN)
    - Radiolysis stability sufficient? => limited stability of DGA
- On-going tests:
  - Use of TK221 (TO-DGA / phosphine oxide) or TK222 (TEH-DGA / phosphine oxide) => TK221 shows higher La retention than DGAs
  - Focus on La/Ac separation
  - Sharp Ac elution => for TK221/2 only in dilute HCl
  - Improved radiolysis stability?
  - Resalting possible? Ac nitrate => Ac chloride
  - Ac data et al: publication under preparation


# Ac-225 separation – method under optimisation

- Two TK221 cartridges for removal of impurities incl. La
- Two additional, optional Pb removal steps (TK102 and TK101)



## Step 1 TK221:


Target dissolved in 2 – 4M HNO<sub>3</sub>  
Ra, Ba, Pb, Sr,... removal with 4M HNO<sub>3</sub>  
Ac elution in ~12M HNO<sub>3</sub> (LNs retained)



## Step 2 TK221:


2x diluted eluate from first TK221  
Rinse with 6M HNO<sub>3</sub> and optional rinses with:  
10M HCl => Bi removal and  
0.05M HNO<sub>3</sub> (Fe, Po removal)  
Ac elution in 0.05M HCl

# TK221 Resin – Ac separation – step one





- In case LN need to be removed
- Two step procedure

- First Ac / LN separation
- First TK221
  - Load from elevated HNO<sub>3</sub>
  - Ac elution in very high HNO<sub>3</sub>
  - LNs, U, Th retained
  - Particular attention to Pb/Sr
    - Elution in 4M HNO<sub>3</sub>




- Second TK221
  - Dilute x2 => load
  - Bi removal 10M HCl
  - Fe removal in 0.05M HNO<sub>3</sub>
  - Ac elution in 0.05M HCl
  - Important: Lanthanides need to be removed upfront (1st TK221)
  - Additional purification on TK101 possible (Ra, Ba, Pb, Sr)

# Optional: TK101 purification step



Optional:  
Ra, Pb, Sr, Ba removal on TK101

Load:  
Ac-225 in <5 mL 0.05M HCl



Ac-225 in 0.05M HCl  
Ra, Pb, Sr remain retained



Data courtesy of B. Russel (NPL)

Optional Pb, Bi, Ra, Sr,... removal step (TK101)  
Pass Ac fraction (0.05M HCl) through TK101  
Ac passes - Ra, Pb, Sr, Bi,... retained




- In case Ra needs to be purified on-column (e.g. dissolved Ra needles) => Use of TK101 for Ra retention / purification
  - Test against Chelex, CEX, TK100
- TK101 => similar to TK100 but ionic liquid replaces HDEHP
  - Both based on same crownether as SR Resin
  - TK100 developed for Sr and Pb uptake also between pH ~2 and 7 (DGT)
    - ⇒ Wagner et al. TK100 discs
    - ⇒ Retains wide range of elements
  - Replacing HDEHP by ionic liquid (=> TK101 Resin) allows for retention of Pb, Sr, Ba, Ra,... from pH ~2 – 7 without extensive extraction of other elements

## TK101 Group 1 and 2





NPL  
National Institute of Standards and Technology





Data provided by  
Russel et al. (NPL)

- Ra retention from water/dilute acid up to  $\sim 0.5$ M  $HNO_3/HCl$
- At higher conc. selectivity closer to SR Resin/TK102 Resin




Data provided by  
Russel et al. (NPL)



- No / extremely low selectivity for TM, Th/U, Ac
- Very strong Pb retention => elution in high HCl or citrate

# Ra separation on TK101




- Good Ra separation when loading from dilute  $\text{HNO}_3/\text{HCl}$
- When eluting Ra in 3M  $\text{HNO}_3$ , Ba, Pb, Sr remain retained
- No retention of U, Th, Pt, Ir,...
- Ra eluted in 3M  $\text{HNO}_3$
- Th and Ba eluted in 8M  $\text{HNO}_3$

# Ra purification


- Work on crown-ether based resin for Ra ongoing
  - Aim: Ra retention from acidic/high  $\text{NO}_3^-$  matrices, high capacity
- In-between: work on TK101 and TK102...

Elution study - 8 mL TK101 Resin




- Ra separation from matrix (e.g. Pt/Ir) and Ba
  - Ra retention and purification at 0.05M  $\text{HNO}_3$
  - Ra elution in 3M  $\text{HNO}_3$
  - Pb, Sr, Ba remain retained
- Might require additional TK102 for Ba removal

Separation on 1 mL TK102 Resin (100 - 200 $\mu\text{m}$ ) - ~0.5BV/min




Separation on 1 mL SR Resin (100 - 150  $\mu\text{m}$ ) - 0.5BV/min



# CU Sheets




- Poster presented at Terachem 2022 (Svedjedet et al.)
- Other than for DGA Sheets not for radionuclidic purity
- QC of Cu radiolabeled peptides (labeled vs free Cu)
  - Shown:  $[^{61}\text{Cu}]\text{Cu-NOTA-octreotide}$
- Spotting/run on three different papers after labeling:
  - Whatman and iTLC without modification and
  - CU extractant impregnated iTLC paper.
- Both iTLC paper (impregnated/non-impregnated) developed in less than 10min, Whatman took 25 – 30 min.
- CU extractant impregnated iTLC paper showed superior resolution
- Beta testing, commercialisation September 2023




- Other systems under development/testing (TK101, ZR,...)


- TO-DGA (normal DGA) and TEH-DGA (branched DGA) impregnated TLC paper
  - Developed at CVUT (Kozempel et al.)
- QC of radionuclides and generator eluents (p.ex. Ra-223, Ac-225/Bi-213, Pb-212, Ge-68/Ga-68 ...)
  - TLC scanner or radiometer/LSC or HPGe after cutting
- Run under acidic conditions => radionuclidic purity



A scheme of chromatographic separation of mixture of  $^{227}\text{Ac}$  and his daughter's nuclides.  $^{227}\text{Th}$  remains on start,  $^{227}\text{Ac}$  has the retention factor ca 0.2,  $^{211}\text{Pb}$  ca 0.7 and  $^{223}\text{Ra}$  ca 0.9.



Radiochromatogram measured immediately after separation. Low abundant radiations of  $^{227}\text{Ac}$  were not detected.



Radiochromatogram measured one hour after separation. Decay and ingrowth of  $^{211}\text{Pb}$  is clearly visible.

- More types of sheets under development (selectivities, geometry, support)
  - ZR, TK201,...
  - 2D TLC for radionuclide screening ?

- Upscale of radiolanthanide separations
- Ra Resins
- Radiometal purification
  - Sb, Pd, Hg, Mn, V, In, Sc, At...
- Improvement of radiolysis stability
- Separation of DTMs
  - Calixarene-based Cs/Rb Resins
  - Mo, Nb, Se, Sb,...
- RP RN in the environment
- In-field preconcentration
  - Impregnated membranes
  - Cartridges
- Passive sampling (DGT)
  - TK100 discs for Sr, Pb, Zn
    - E.g. [Wagner et al.](#): Labile Pb and Sr in soil samples via DGT
  - CL resin for iodine, CA for Ra,...
- Rapid tests
  - Range of impregnated PSm resins
    - Uni Barcelona
  - Impregnated membranes
  - Range of 'Test sticks'
    - NPL, JCU
- Other 'geometries' & 'Non-resin' separation materials
- Microfluidics





Thank you for your attention!



**SUBSCRIBE TO OUR NEWSLETTER**

To keep updated with our latest developments, news and agenda for a year, subscribe to the TrisKem Info here

