

Development of a Sodium Bismuthate-Coated Polyacrylonitrile Resin for the Separation of Oxidized Americium from Curium

Samantha A. Labb
65th RRMCM – TrisKem Workshop | November 3, 2022




Nuclear Data Needs for Stockpile Stewardship

- NNSA 2021 Stockpile Stewardship and Management Plan
 - Ensure safety and reliability of the stockpile without underground testing
- Nuclear Data Needs
 - Accurate cross section measurements
 - Am-242m/Am-242g Isomeric Ratio





Measurement Strategy



Required: Complete separation and recovery of Am and Cm


Am/Cm SEPARATION USING 5% NaBiO₃ IN CELITE



NaBiO₃ OXIDIZING AGENT

Ilmenite Structure → Ion Exchange Properties


- Octahedral bismuth centers and sodium cations
- Literature shown for alkali, alkaline earth, and U(VI)



NaBiO₃-PAN RESIN DEVELOPED BY TrisKem

- Resin Preparation

- NaBiO₃ not bound to resin
- Not prepared with diluents like EXC resins




- Particle Size Analysis

- Dried resin, broke up clumps, and separated with hand sieve
- 125 - 500 micron range




Effect on $[\text{HNO}_3]$

- Surface of sodium bismuthate is basic
- Change in Acid Concentration
 - HNO_3 solutions titrated prior to contact with resin
 - 1.5 mL of acid added to 50 mg of resin x2
 - Post-contact titrations
- Acidification of surface necessary




Hydrolysis Study

1. Contacted 1.5 mL of 0.01 and 0.05 M HNO_3 with 50 mg NaBiO_3 for 1 hour
2. Separated by centrifugation, transfer 1.45 mL supernatant to empty microcentrifuge tube
3. Spike with 50 μL radionuclide stock
4. Shake for 1 hour
5. Analyze 1 mL aliquot via LSC



METHOD: BATCH CONTACT STUDIES



1. Preconditioning

- 50 mg resin
- Desired acid concentration
- Agitate for 1 hour


2. Spiking

- $^{241}\text{Am}/^{244}\text{Cm}$
- Agitate for 1 hour


3. Analysis

- Filter sample using PTFE syringe filter
- Transfer aliquot of the filtered eluent for analysis


Am/Cm on NaBiO₃ Solid (Idaho National Lab)




75 wt% NaBiO₃-PAN Concentration Dependency




75 wt% NaBiO₃-PAN Precondition Time Dependency (0.01 M HNO₃)




75 wt% NaBiO₃-PAN Contact Time Dependency (0.01 M HNO₃)




75 wt% NaBiO₃-PAN Temperature Dependency (0.01 M HNO₃)



Am/Cm on PAN Beads




METHOD: CHROMATOGRAPHIC SEPARATIONS



COLUMN OBSERVATIONS



Small floaters had to be decanted when prepping slurry to prevent column clogging



Resin turning black when acid introduced ($\text{Bi(V)} \rightarrow \text{Bi(III)}$)
Kozma, *et. al. Solid State Chem.* **2018**, 263, 216-223)




Gas production restricting flow rate under gravity




Color change 1 week post column

ATTEMPT #1: NaBiO₃-PAN COLUMN




- 2% Am-241 recovery for every 10 mL of 0.01 M HNO₃ eluted
- 100% Am and Cm recovery in 0.05 M HNO₃
- Retention is **TOO HIGH** at 0.01 M HNO₃
- Very restricted flow under gravity, gas production


COLUMN OBSERVATIONS



Density differences and wide range of particle sizes affect flow dynamics and reproducibility




Better distribution with the 10 wt% NaBiO₃-PAN column. Involved a lot of mixing with pipette.




Dry packing shows more promising bed uniformity – tight packing released remaining solution in resin (pH ~11)

ATTEMPT #2: 5% NaBiO₃-PAN IN PAN-Normal




5 wt% NaBiO₃-PAN
95 wt% PAN-Normal




- Non-homogenous mixing of the two resins due to density
- Uneven bed packing
- 98% Am recovery, 100% Cm recovery
- Retention is **TOO LOW** at 0.01 M HNO₃




Attempt #3: 10% NaBiO₃-PAN with PAN Beads



10 wt% NaBiO₃-PAN
90 wt% PAN-Normal




- 90% Am and 71% Cm recovery
- Separation! ... with some Am breakthrough in the Cm fraction
- Starting to find the sweet NaBiO₃ concentration spot




Next... Improving the resin

- At high amounts of NaBiO_3 coating, both radionuclides are strongly retained
 - Vary NaBiO_3 loading on resin
 - 10%, 25%, 50%, 75%
- Exploring PES beads as backbone
- Comparing NaBiO_3 sources & purity


Am/Cm on 75% NaBiO₃-PAN (New Lot)




Am/Cm 75 wt% NaBiO₃-PAN (Old Lot)




Am/Cm 75 wt% NaBiO₃-PAN (New Lot)



NaBiO₃ Solid (Idaho National Lab)




NaBiO₃ Solid (TrisKem - As Received)




NaBiO₃ Source Purity - Peroxides

- NaBiO₃ Synthesis
- Peroxide contaminants → incomplete or non-oxidation of Am
 - Am(III) → Am(VI) → Am(V) + Am(III)
 - Am(III) → Am(V)
- Remove through repeated washes with DI-H₂O (?)


NaBiO₃ Solid (TrisKem - As Received)




NaBiO₃ Solid (TrisKem - Washed)




Am/Cm on 50 wt% NaBiO₃-PAN



Am/Cm on 25 wt% NaBiO₃-PAN



Am/Cm on 10 wt% NaBiO₃-PAN




Polyacrylonitrile Support


- No evidence of oxidation
 - Identical Am/Cm behavior → trivalent speciation
 - Organic PAN reducing Am(V)/Am(VI)?
- No evidence of ion exchange
 - No Cm(III) adsorption → exchange site access
 - NaBiO₃ on surface or in pores?
 - PAN incorporation into solid structure?
 - Zirconium Phosphate

Polyethersulfone (PES) Support

- 10 wt% NaBiO₃ coating
- Greater stability in acid than PAN
- Potential benefit of additional electrostatic interactions



Am/Cm on 10% NaBiO₃-PES Fresh



Alternative Chromatographic Systems

- Inorganic Supports/Filter Aids
 - Silica Gel/Powder → Mixed Bed or Coated
 - Alumina
- Pre-treatment with NaBiO_3 → Cation Exchange
 - Zirconium Phosphate
 - Zirconium Oxide

NaBiO₃ Dissolution/Degradation

- What is the concentration of total Bi (Bi(V) + Bi(III)) in solution as a function of:
 - [HNO₃]
 - Time
 - Radionuclide concentration, size, charge
- ICP-OES: Total [Bi]
- UV-Vis-NIR: Bi(V):Bi(III)
- **Bi(III) Interference: How is Cm(III) adsorption affected?**

Beyond Resin Development

- Oxidation & Speciation
 - UV-Vis analysis in the presence of NaBiO_3 to determine speciation and oxidation kinetics
 - Function of time, $[\text{NaBiO}_3]$
 - Compare speciation determination methods – fluoride precipitation vs. UV-Vis
- Ion Exchange
 - What is the ion exchange mechanism for Am and Cm with Na^+ ?
 - Potential complexation vs ion exchange??
 - EXAFS / XANES Studies

Acknowledgements

- Dr. Ralf Sudowe
- Dr. Tom Johnson
- Dr. Alex Brandl
- Dr. Chuck Henry
- Aude Bombard, TrisKem Int.
- Steffen Happel, TrisKem Int.
- Dr. Evelyn Bond
- Dr. Todd Bredeweg
- ERHS
- G. T. Seaborg Institute



**MOUNTAIN
& PLAINS ERC**

*A NIOSH Education and Research Center for
Occupational & Environmental Health & Safety*

This research was supported by the grant T42OH009229, funded by the National Institute of Occupational Safety and Health in the Centers for Disease Control and Prevention. Its contents are solely the responsibility of the authors and do not necessarily represent the official views of the Centers for Disease Control and Prevention or the Department of Health and Human Services.



This research was supported under award NRC-HQ-60-17-G-0005 from the U.S. Nuclear Regulatory Commission. The statements, findings, conclusions, and recommendations are those of the author(s) and do not necessarily reflect the view of the US Nuclear Regulatory Commission.